Skill Sheet

Nasopharyngeal Airway

Objective: DEMONSTRATE the recommended procedure for insertion of a Nasopharyngeal Airway.

References:

PHTLS (Military Version) Eighth Edition, Jones and Bartlett Learning ATLS Manual, American College of Surgeons, Chicago, IL

Evaluation: Students will be evaluated as a Pass/Fail (P/F). The instructor will verify the accuracy of the student's ability to insert a Nasopharyngeal Airway on an airway trainer by means of observing the student's procedures and technique.

Materials:

Student Checklist

Airway Simulator, Nasopharyngeal Airway, surgical lubricant

Instructor Guidelines:

- 1. Provide each instructor with a Student Checklist.
- 2. Ensure student has all student-required materials.
- 3. Read the Learning Objective and the evaluation method to the student.
- 4. Explain the grading of the exercise.
- 5. Allow time for the students to extract the information required from the instructor-provided scenario.

Performance Steps:

- 1. Assemble and test all necessary equipment.
- 2. Assess the upper airway for visible obstruction.
- 3. Open the airway with a chin lift/jaw thrust maneuver.
- 4. Verbalize the indications for a Nasopharyngeal Airway unconscious patient.
- 5. Lubricate the Nasopharyngeal Airway with a surgical lubricant.
- 6. Insert the airway into the nose at a 90-degree angle to the face. Avoid aiming upwards towards the top of the head. Insert all the way to the flange.
- 7. Use a rotary and/or back-and-forth motion to facilitate insertion.
- 8. If unable to insert on one side of the nasal passage, take it out and try the other side.

Insert a Nasopharyngeal Airway

Completed

Task

Evaluator:

1 ask		Jumpiete	u
	1st	2nd	3rd
Assessed the upper airway for visible obstruction.	P / F	P / F	P / F
Opened the airway with a chin lift/jaw thrust maneuver.	P / F	P / F	P / F
Lubricated the Nasopharyngeal Airway with a surgical lubricant.	P / F	P / F	P / F
Inserted the airway into the nose at a 90-degree angle to the face. Avoided aiming upwards towards the top of the head. Inserted all the way to the flange	P / F	P / F	P / F
Used a rotary and/or back-and-forth motion to facilitate insertion.	P / F	P / F	P / F
If unable to insert on one side of the nasal passage, took the NPA out and inserted on the other side.	P / F	P / F	P / F
ritical Criteria: Aimed the NPA towards the top of the head during insertion Performed procedure in a manner that was dangerous to the			
nator's Comments:			

Student Name: _____ Date: _____

Pass: ____ Fail: ____